

# Prodigy™ SCS System



The Prodigy MRI™ Implantable Pulse Generator (IPG) is a rechargeable MR Conditional\*\* device featuring both BurstDR™ stimulation and tonic stimulation modes, providing more patients with the opportunity for maximum therapy effectiveness. Now patients who may require access to an MRI scan\*\* can also access the benefits of BurstDR stimulation—a proven, patient preferred option to manage chronic pain.<sup>1,6</sup>

## **BURSTDR STIMULATION: A PROVEN TECHNOLOGY**

Supported by a body of clinical evidence, BurstDR stimulation has been proven to deliver superior pain relief over tonic stimulation.<sup>1,6</sup> It works by emulating natural firing patterns in the brain<sup>7</sup> and is believed to modulate both the sensory and emotional pathways—giving patients relief from both their pain sensation and the suffering\*\*\* associated with the pain.<sup>1,6,7</sup>

## **SMALLEST DEVICE**

The Prodigy MRI IPG is the smallest rechargeable IPG available.<sup>1-5</sup> Smaller devices may reduce wound closure time and improve cosmesis in patients.

## **LONG BATTERY LIFE**

The device features the longest projected battery life,<sup>1-5</sup> which includes 10 years of practical recharging.<sup>1</sup>

## **UPGRADEABLE THERAPY**

The Prodigy MRI IPG allows for upgradable technology innovations upon approval—without the need for surgical revision. The Prodigy MRI IPG is MR Conditional for Head and Extremity scans. For scan details of our current neuromodulation devices, refer to the [resources for radiology professionals](#).

## **MOVE MORE PATIENTS FROM PAIN TO RELIEF**

The Prodigy MRI™ IPG is part of our spinal cord stimulation (SCS) portfolio. Our vision is to transform the management of chronic pain by offering advanced technologies

designed to improve patient outcomes and experience to bring the benefits of SCS to more patients. Read more about the [portfolio](#) and our approach to [chronic pain management](#).

Table 1. IPG specifications

Model	3772	3799
MRI status	MR conditional	Untested
Height	4.8 cm (1.89 in)	
Length	5.3 cm (2.09 in)	
Thickness	0.95 to 1.1 cm (0.37 to 0.43 in)	
Weight	29.0 g (1.0 oz)	
Volume	17.7 cm <sup>3</sup> (1.08 in <sup>3</sup> )	
Power source	Rechargeable lithium ion cell	
Storage temperature	-10°C–55°C (14°F–131°F)	
Storage humidity	10%–90% (noncondensing)	
Storage pressure	70–150 kPa (10.2–21.8 psi)	
Connector strength	Exceeds EN45502-1 requirements	




Table 2. Operating parameters for the IPG

Parameter	Tonic Range	Tonic Steps	Burst Range*	Burst Steps*
Pulse width	50–500 $\mu$ s	Alternating 12 and 13 $\mu$ s (starting with 12 $\mu$ s)	50–1000 $\mu$ s	50 $\mu$ s
Frequency	2–200 Hz	2 Hz	—	—
	200–500 Hz	10 Hz	—	—
	500–1200 Hz	20 Hz	—	—
Burst rate frequency	—	—	10–60 Hz	10 Hz
Intraburst frequency	—	—	250–500 Hz	10 Hz
	—	—	500–1000 Hz	20 Hz
Amplitude	0–25.5 mA (max 12 V)	0.1–1.0 mA	0–12.75 mA	0.05–0.50 mA

NOTE: Columns with \* represent operating parameters for BurstDR™ stimulation programs on IPGs capable of BurstDR stimulation mode.

NOTE: The number of stim sets in use for a tonic program governs the maximum frequency (1200/number of stim sets).

NOTE: The maximum current depends on the impedance, frequency, and pulse width settings.

## REFERENCES

\*BurstDR™ stimulation, patented technology exclusively from St. Jude Medical, is also referred to as Burst stimulation in clinical literature.

\*\*Within approved scan parameters. Refer to the Instructions for Use for full details of the St. Jude Medical™ Prodigy MRI™ IPG MR Conditional scan parameters.

\*\*\*Pain and suffering as measured by VAS.

1. St. Jude Medical. (2016). St. Jude Medical™ Prodigy™ Neurostimulation System Programming and Reference Manual. Plano, TX.
2. Medtronic. (2011). RestoreSensor™. Multi-program Rechargeable Neurostimulator Implant Manual. Minneapolis, MN.
3. Boston Scientific. (2011). Precision™ Spinal Cord Stimulation System Clinician Manual. Valencia, CA.
4. Nevro Corporation. (2012). Nevro Physician Implant Manual 10186-Eng Rev. F. Menlo Park, CA.
5. Nuvectra. (2014). Algovita™. Spinal Cord Stimulation Patient System Manual. Plano. TX.
6. St. Jude Medical. (2016). St. Jude Medical™ Proclaim™ Neurostimulation System Clinician's Manual. Plano, TX.
7. De Ridder, D., Vanneste, S., Plazier, M., & Vancamp, T. (2015). Mimicking the brain: Evaluation of St. Jude Medical's Prodigy Chronic Pain System with Burst Technology. *Expert Review of Medical Devices*, 12(2), 143–150. <http://dx.doi.org/10.1586/17434440.2015.985652>